

TOUCH 5

Electro Surgical Generator



Quality, Performance & Reliability!

A **NEW** GENERATION OF ELECTROSURGICAL GENERATORS

TOUCH 5 SERIES OF ELECTROSURGICAL GENERATORS WITH:

- 5-inch LCD touch screen display
- Improved vessel-sealing technology
- Smoother and faster coagulation
- Superior spray coagulation
- First-in-the-segment Bipolar in Saline (Bi-TUR) variant



IMPROVED VESSEL SEALING

- Enhanced and improvised vessel sealing performance with BISSECTOR
- Faster vessel sealing, lower vapour generation and lower lateral thermal spread than ELSY-360 M+

BIPOLAR IN SALINE VARIANT

- Experience ALAN's exceptional Plasma Bipolar Energy in a more economical package

ILLUMINATED SOCKETS

- Illuminated connector sockets for better visibility in challenging light conditions

WARRANTY

- First-in-the-segment warranty terms of 2-years for all Touch5 variants

CAPACITIVE RETURN ELECTRODE

- Improved support for Capacitive Return Electrode, for unparalleled safety of patients
- Gold standard for patient safety with respect to prevention of patient plate site burns

PPCM & TISSUE MONITORING TECHNOLOGY

- Patient Plate Contact Monitoring technology for enhanced safety
- Continuous monitoring of tissue impedance for dynamic control of current and voltage

MULTIPLE NEW VARIANTS OF TOUCH 5



COMBI	COMBI-PLUS	ALLIGATURE	PLASMA-PLUS	PLASMA
MONOPOLAR	MONOPOLAR	MONOPOLAR	MONOPOLAR	
BIPOLAR	BIPOLAR	BIPOLAR	BIPOLAR	BIPOLAR
		VESSEL SEALING		
			BI-TUR	BI-TUR

CAPACITIVE RETURN ELECTRODE



- Gold standard for patient safety w.r.t. prevention of return electrode site burns
- No hair removal, no skin tears, and no adhesive allergic reactions
- Connection ports available on all 4 corners for easy access from electro surgical generator on all sides in OR

ALAN

Manufactured by:

ALAN ELECTRONIC SYSTEMS PVT. LTD.

Factory & Head Office:

Plot No. F-75, MIDC, Additional Ambarnath Industrial Area, Anand Nagar, Jambivali, Ambarnath (East) - 421506, District : Thane, Maharashtra, INDIA. Tel : +91-251-2628000.

Website: www.alanindia.net

India:

+91 93205 53969

sales.medical@alanindia.net

International:

+91 88790 03067

export@alanindia.net

