




Product feature

- Air temp servo-controlled by computer
- Various and self-check alarms
- Removable humidity reservoir, easy to clean
- With the transfusion shelf and tray
- RS-232 connector
- Phototherapy unit (option)



Tube phototherapy unit(option)



Monitor tray



The humidity reservoir is removable



LED Phototherapy unit(option)


Specifications Infant Incubator

| Model                      | BIN-3000A(G)               |
|----------------------------|----------------------------|
| Electrical requirement     | ~ 220V 50Hz or ~ 230V 60Hz |
| Power consumption          | ≤ 450VA                    |
| Air mode temp. range       | 25.0℃ ~ 37.0℃              |
| Skin mode temp. range      | -                          |
| Air over temp. range       | < 38.0℃                    |
| Air mode high temp. range  | -                          |
| Skin mode high temp. range | -                          |
| High temp. mode alarm      | -                          |
| Air mode alarm             | ± 3.0℃                     |
| Skin mode alarm            | -                          |
| Skin sensor accuracy       | -                          |
| Air control accuracy       | ≤ 0.5℃                     |
| Mattress temp. uniformity  | ≤ 0.8℃                     |
| Noise                      | ≤ 55dB(A)                  |



Product feature

- Air temp. servo-controlled by computer
- > 37℃ temperature setting
- Various and self-check alarms
- Removable humidity reservoir, easy to clean
- With the transfusion shelf and tray
- RS-232 connector
- Phototherapy unit (option)



Phototherapy unit(option)  
BBP-1000B/2000B



Monitor tray



The humidity reservoir is removable



Phototherapy unit(option)  
BBP-3000B

Specifications Infant Incubator

| Model                     | BIN-3000A                  |
|---------------------------|----------------------------|
| Electrical requirement    | ~ 220V 50Hz or ~ 230V 60Hz |
| Power consumption         | ≤ 450VA                    |
| Air mode temp. range      | 25.0℃ ~ 37.0℃              |
| Air over temp. range      | < 38.0℃                    |
| Air mode high temp. range | 37.1℃ ~ 39.0℃              |
| High temp. mode alarm     | < 40.0℃                    |
| Air mode alarm            | ± 3.0℃                     |
| Air control accuracy      | ≤ 0.5℃                     |
| Mattress temp. uniformity | ≤ 0.8℃                     |
| Noise                     | ≤ 55dB(A)                  |