

BRW-4000A/4000B



| Model                   | BRW-4000A | BRW-4000B |
|-------------------------|-----------|-----------|
| screen                  | LCD       | LCD       |
| skin temperature sensor | ●         | ●         |
| suction                 | —         | ●         |
| O <sub>2</sub> bottle   |           | option    |
| electronic lift         |           | option    |
| phototherapy            |           | option    |

being

BRW-4000B



LCD Touch Screen

- Shows all data including temperatures, APGAR timer, time of phototherapy
- 72 hours data record
- Various alarms



Infrared Ceramic Radiator

- Explosion-proof, safety assurance

LED Phototherapy Unit



Mattress Tilting Adjust

- Free-step adjustment  $\pm 12^\circ$



X-ray Tray under the bed



Suction

Specifications

| Model                     | BRW-4000A                    | BRW-4000B |
|---------------------------|------------------------------|-----------|
| Electrical requirement    | ~ 220V 50Hz or ~ 230V 60Hz   |           |
| Power consumption         | $\leq 700\text{VA}$          |           |
| Air mode temp. range      | 25.0℃ ~ 37.0℃                |           |
| Over temp. alarm          | $< 38.0^\circ\text{C}$       |           |
| High mode temp. range     | 37.1℃ ~ 38.0℃                |           |
| High mode temp. alarm     | $< 39.0^\circ\text{C}$       |           |
| Temp. deviation alarm     | $\pm 1.0^\circ\text{C}$      |           |
| Skin sensor accuracy      | $\leq \pm 0.3^\circ\text{C}$ |           |
| Mattress temp. uniformity | $\leq 2.0^\circ\text{C}$     |           |
| Heating rate              | $\leq 30\text{min}$          |           |
| APGAR timing              | 1min, 5min, 10min            |           |