

BIOBASE BPR-5V298

- Laboratory Refrigerator (2~8°C)
- Glass door version.
- This is single tall door with 5 shelves.

SPECIAL FEATURES FOR BPR-5V298

1. For Lab / Pharmacy / Vaccine Storage:

2. High Visibility, No Fog:

- Internal LED lighting auto on/off with door
- You can read stock without opening door = temperature stays stable

3. Accurate Forced Air Cooling:

- Air-cooled fin
- Auto defrost
- Rapid recovery after door opening

4. Intelligent Alarm System (Audible + Visual):

- High temperature alarm / Low temperature alarm
- Power failure alarm - with backup battery
- Sensor failure / Thermostat failure alarm

5. Data & Security for Lab Accreditation:

- Microprocessor control with Large LED digital display
- USB port for temperature data download
- Door lock with key

6. Construction:

- Internal: Embossed aluminum sheet easy to clean, anti-corrosion
- External: Cold-rolled steel with anti-bacteria coating
- 4 universal casters with brakes