


BJS - 004

C-ARM COMPATIBLE SEMI ELECTRIC OT TABLE



Salient Features

- Electric Up-down by Cable Remote • Eccentric Pillars • Stainless Fitting • Detachable, Head, Leg & Pelvic Section • Four Sectional Radio-Translucent Top • Base & Cover Made of Attractive Easy to clean Stainless Steel • Sophisticated Mechanics Provide Smooth, Step less articulation of the table top for precise patient positioning. • Use Remote for Height Adjustment • In Built Kidney Bridge.

Suitable For

General Surgery
Urology Surgery
Gynecology
Laparoscopy

Vascular
Cardiac Surgery
Neurology Surgery
Plastic Surgery

Standard Accessories

Anesthetic Screen
Shoulder Supports with Pad
Arm Boards with Pad
Lateral Support with Pad
Knee Crutches
Water Proof Mattress

Technical Data

Length of Table	1910 mm
Width of Table Top	535 mm
Height (Min- Max)	750 - 1050 mm
Trendelenburg (Rev.)	25°
Lateral Tilt	±20°
Back Section	±15° - 90°
Head Section	±45° - 90°
Leg Section	±45° - 90°
Floor Locking System	Manual
Battery Backup	Optional