


C-ARM COMPATIBLE 6 FUNCTION FULL ELECTRIC OT TABLE



Salient Features

THIS FULL ELECTRIC TOP SLIDE. O.T Table is suitable for all kind of surgeries i.e. general surgery, cascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others. The base & column of the table is made of high quality medical stainless steel (Grade - 304). Non-reflecting surface, antibacterial & easy to clean. We provide with a top sliding of 5 inch back/forward. Height adjustment, Lateral Tilts, Trendelenburg, and Flex - Refl are realize by Remote control.

Suitable For

General Surgery
Urology Surgery
Gynecology
Laparoscopy
Vascular
Cardiac Surgery
Neurology Surgery
Plastic Surgery

Standard Accessories

Anesthetic Screen
Shoulder Supports with Pad
Arm Boards with Pad
Lateral Support with Pad
Knee Crutches
Water Proof Mattress

Steel Quality
304 Medical Grade Stainless Steel

Technical Data

Length of Table	1910 mm
Width of Table Top	535 mm
Height (Min- Max)	750 - 1000 mm
Trendelenburg (Rev.)	25°
Lateral Tilt	±20°
Back Section	±15° - 90°
Head Section	±45° - 90°
Leg Section	±45° - 90°
Floor Locking System	Remote
Top Slide	8" Reverse & 8" Forward

 ADDITIONAL CONTROL PANEL

 MEMORY FUNCTION

 ZERO POSITION

 AUTO FLOOR LOCKING