

BJS - 007

ELECTRIC OPHTHALMIC OT TABLE



Length
1900 mm



Height
750 - 1000 mm



Width of Table
520 mm



Trendelenburg
25°



Salient Features

This Ophthalmic Table is specially designed for ophthalmic surgery. Three sectional table is provided with drainage system in the head section. The table has sturdy base, column made up of stainless steel, antibacterial and is easy to clean. Electromotive height and Trendelenburg positions with the help of light weight electric foot control and with remote control.

Standard Accessories

- Anesthetic Screen
- IV Pole
- Hand Rest Mattress

Technical Data

Length of Table	1910 mm
Width of Table Top	535 mm
Height (Min- Max)	750 - 1000 mm
Trendelenburg (Rev.)	25°
Trendelenburg	25°
Leg Section	0° - 90°

Steel Quality
304 Medical Grade Stainless Steel