


BJS - 008

GYNECOLOGICAL OBSTETRIC TABLE

OT TABLE



Salient Features

Obstetrics operating room table perfect for all kinds of gynecological operations. The cover for base and column is made of 304 stainless steel, good performance on anti-rust and easy to clean THIS FULL ELECTRIC TABLE WITH OPTIONAL BATTERY BACKUP OT. OT Table is suitable for all kinds of surgeries i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others. The base & column of the table is made of high quality medical stainless steel (Grade - 304). Non-reflecting surface, antibacterial & easy to clean. We provide in this table Lateral Tilts, Trendelenburg, and Flex - Reex, Chair Position and Longitudinal are realize by Remote control.

Suitable For

General Surgery	Vascular
Urology Surgery	Cardiac Surgery
Gynecology	Neurology Surgery
Laparoscopy	Plastic Surgery

Steel Quality
304 Medical Grade Stainless Steel

Standard Accessories

Birth Handle 2pcs
Gynae Tray
Anesthetic Screen
Shoulder Supports with Pad
Arm Boards with Pad
Lateral Support with Pad
Knee Crutches
Water Proof Mattress

Technical Data

Length of Table	1910 mm
Width of Table Top	535 mm
Height (Min- Max)	750 - 1000 mm
Trendelenburg (Rev.)	25°
Lateral Tilt	±20°
Back Section	±15° - 90°
Head Section	±45° - 90°
Leg Section	±45° - 90°
Floor Locking System	Manual
Battery Backup	Optional