




BBP-2000B

BBP-2000B(LED)

- Easily move with wheels
- LCD display shows time information, lighting intensity
- With white LED lights for easy examination
- Programmable system controller to set different lighting intensity for different time


BBP-3000B(LED)



BBP-3000B

Specifications

Model	BBP-2000B	BBP-3000B
Bilirubin radiancy (40cm) Bilirubin radiancy (50cm)	3750 μ w/cm ² (62.5 μ w/cm ² /nm) 3200 μ w/cm ² (53 μ w/cm ² /nm)	
Lighting intensity	33%,67%,100%	
White light	Yes	
Power	≤ 60VA	
Wave length	400nm ~ 550nm	
Radiant space	400 × 250mm	
Radiant Uniformity	> 0.4	
Therapy time set	0~99h59m	
Height range	1350mm ~ 1650mm	—



**BBP-400A(4tubs)
BBP-500A(5tubs)**

- Available stand pole up-down adjustment
- Available angle adjustment 0° , 30° , 60°
- Wave length: 400-500 um
- Fitted with any infant bed and incubator
- 4 tubes or 5 tubes

BBP-1000B(LED)

Phototherapy unit can be put on the top of incubator



BBP-1000B

Specifications

Model	BBP-1000B
Bilirubin radiancy (40cm) Bilirubin radiancy (50cm)	3750 μ w/cm ² (62.5 μ w/cm ² /nm) 3200 μ w/cm ² (53 μ w/cm ² /nm)
Lighting intensity	33%,67%,100%
White light	Yes
Power	≤ 60VA
Wave length	400nm ~ 550nm
Radiant space	400 × 250mm
Radiant Uniformity	> 0.4
Therapy time set	0~99h59m